

Achieving Optimal Connectivity in the Histology Lab: Using LANs and WANs to Improve Quality and Productivity

Anthony St-Onge, MLT, Lean Six Sigma Green Belt Strategic Business Manager, West Region

Christopher Marino, MBA, Lean Six Sigma Green Belt Strategic Business Manager, East Region



Agenda

- Connectivity Overview
- Customer Example
- Suggestions for Implementation
- Q & A / Discussion











<u>Criteria</u>

- Quality
- Patient Safety
- Improve work process
- Ease of use
- Flexibility
- Cost

Challenges

- Increasing workload
- Improving TAT
- Maximizing use of resources
- Improve outcomes
- Traceability
- Space limitations



Why connectivity?



- Eliminate manual data entry
- Reduce possibility of errors
- Maximize employee resources
- Track specimens and slides
- Manage workflow / workload
- Growth without increase in FTEs
- Better patient care!



Behind every stain there is a face





What is connectivity?



"A Program or Device's ability to <u>link</u> with other Programs and Devices for <u>exchange</u> (Import / Export) of data in various formats"





What is LAN and WAN connectivity?



LAN - Local Area Network

A LAN connects network devices over a relatively short distance. A networked office building, school, or home usually contains a single LAN, though sometimes one building will contain a few small LANs (perhaps one per room), and occasionally a LAN will span a group of nearby buildings.

WAN - Wide Area Network

As the term implies, a WAN spans a large physical distance and is a geographically dispersed collection of LANs. A network device called a router connects LANs to a WAN.



Now available for histology... ... connectivity options to meet your needs





Single server with LAN and/or WAN



Single Server



Critical information can be accessed where it is most needed, when it is most needed.



Possibilities for histology





- LAN seat deployment enables info access from any computer
- Print on demand labels from LIS or individual workstations
- Standardization of staining sharing of protocols across instruments
- Seamless menu and instrument expansion and modifications
- Centralized antibody preparation, optimizations, and validations
- Inter-lab QC proficiency management
- Multi-site contingency for instruments and reagents
- Track and trace slide status, lot #s, reagent expiry
- Reporting functions easy data mining, monitor usage, control material
- Share reagents, send pre-labeled slides to other sites

WAN connectivity





A Better Path[™]

Other possibilities for histology





Remote customer support

- Perform software upgrades as soon as they are available
- Troubleshoot problems without travel delays
- Assist with data collection



Track and trace solutions on LAN





Know where every specimen, cassette, block, or slide is in the process from the time it is received until the slide is archived

WAN for multi-site setting





Monitor remotely or from anywhere within the healthcare system

Every organization is unique













- Population = 1.2 million
- Size = 649,950 km² (250,950 sq mi)



Utilized connectivity solutions at DSM













DSM IHC instrumentation





LAN and WAN connectivity benefits





DIAGNOSTIC SERVICES MANITOBA

- Reduction of potential for errors = Quality
- Optimized human resources utilization
- Save technical time
- Pathologist satisfaction
 - Flexibility to adapt menu
 - Standard staining
- Cost reduction
- Improved TAT





Wide Area Network (WAN) utilized applications at:



DIAGNOSTIC SERVICES MANITOBA

LAN seat utilization







- Agilent Technologics











Wide Area Network (WAN) utilized applications at:



DIAGNOSTIC SERVICES MANITOBA

- LAN seat utilization
- Print on demand at Grace Hospital





Print on demand at Grace Hospital







Wide Area Network (WAN) utilized applications at:



DIAGNOSTIC SERVICES MANITOBA

- LAN seat utilization
- Print on demand at Grace Hospital
- Centralized antibody preparation, optimization, and validations



Centralized preparation, optimization and validations







Wide Area Network (WAN) utilized applications at:



DIAGNOSTIC SERVICES

- LAN seat utilization
- Print on demand at Grace Hospital
- Centralized antibody preparation, optimization, and validations
- Seamless menu expansion and modifications



Seamless menu expansion







Wide Area Network (WAN) utilized applications at:



DIAGNOSTIC SERVICES

- LAN seat utilization
- Print on demand at Grace Hospital
- Centralized antibody preparation, optimization, and validations
- Seamless menu expansion and modifications
- Standardization of staining



Standardization of staining







Wide Area Network (WAN) utilized applications at:



DIAGNOSTIC SERVICES

- LAN seat utilization
- Print on demand at Grace Hospital
- Centralized antibody preparation, optimization, and validations
- Seamless menu expansion and modifications
- Standardization of staining
- Multi-site contingency: Instruments, reagents

Multi-site contingency: Instrumentation





Multi-site contingency: Reagents







Wide Area Network (WAN) utilized applications at:



DIAGNOSTIC SERVICES

- LAN seat utilization
- Print on demand at Grace Hospital
- Centralized antibody preparation, optimization, and validations
- Seamless menu expansion and modifications
- Standardization of staining
- Multi-site contingency: Instruments, reagents
- Track and trace: Slides status



Track and trace: Status of slides









Wide Area Network (WAN) utilized applications at:



DIAGNOSTIC SERVICES

- LAN seat utilization
- Print on demand at Grace Hospital
- Centralized antibody preparation, optimization, and validations
- Seamless menu expansion and modifications
- Standardization of staining
- Multi-site contingency: Instruments, reagents
- Track and trace: Slides status
- Reports: Easy data mining, monitor usage





Protocol steps Slid	le log	
Time		Log
11/21/2012 7:15:46	5 AM	PT Link run started on PT 2 at 11/21/2012 7:15 AM
11/21/2012 7:56:14	1 AM	PT Link run completed on PT 2 11/21/2012 7:56 AM.
11/21/2012 7:56:14	1 AM	TR Fluid was High LHSC mixed 11/20/2012 Cycle 1 of 3
11/21/2012 7:56:14	1 AM	PT Link run log: "2012-11-21-0715-PT 2-L.dta"
11/21/2012 8:16:09	AM 6	Staining started
11/21/2012 8:17:21	AM	Rinsed with buffer (00085374 2012-12-28)
11/21/2012 8:24:33	3 AM	Blow and apply 150 μl EnVision FLEX Peroxidase-Blocking Reagent (00085345 2013-08-31) [1,3]
11/21/2012 8:29:40) AM	Rinsed with buffer (00085374 2012-12-28)
11/21/2012 8:48:34	1 AM	Blow and apply 150 µl FLEX Cytokeratin, MxH (AE1/AE3) (10067833 2014-02-28) [1,3]
11/21/2012 9:08:40) AM	Rinsed with buffer (00085374 2012-12-28)
11/21/2012 9:38:40) AM	Rinsed with buffer (00085374 2012-12-28)
11/21/2012 9:45:07	7 AM	Blow and apply 150 µl EnVision FLEX /HRP (00085348 2013-08-31) [1,3]
11/21/2012 10:05:1	7 AM	Rinsed with buffer (00085374 2012-12-28)
11/21/2012 10:30:2	20 AM	Rinsed with buffer (00085374 2012-12-28)
11/21/2012 10:39:0	01 AM	Blow and apply 150 µl Substrate Working Solution (mix) (00086495 2013-08-31) [1,3]
11/21/2012 10:49:0)2 AM	Rinsed with buffer (00085374 2012-12-28)
11/21/2012 10:49:0)2.AM	Staining completed



Reagent Bottle Information



Usage

Date	Amount in µl
11/19/2012	19800

Slides this reagent was applied to

	Case number	Name
107383	S12-45954	
107384	S12-45954	
107385	S12-45954	
107386	S12-45954	
107387	S12-45954	
107388	S12-45954	
107706	S12-45343	
107707	S12-45343	
107708	S12-45343	
107709	S12-45343	
107710	S12-45343	
107711	S12-45343	
107712	S12-45343	
107719	S12-45824	
107720	S12-45824	
107721	S12-45824	
107722	S12-45824	
107723	S12-45824	_
107724	S12-45824	

Protocol	Completion time
Toxoplasma	11/19/2012 11:23:58 AM
CD20	11/19/2012 11:24:04 AM
CD79	11/19/2012 11:24:11 AM
BCL2	11/19/2012 11:24:18 AM
UNC Mouse	11/19/2012 11:24:25 AM
UNC Rabbit	11/19/2012 11:25:17 AM
CKAE1/AE3	11/19/2012 1:38:34 PM
CKCAM5.2	11/19/2012 1:38:41 PM
SMA	11/19/2012 1:39:33 PM
Desmin	11/19/2012 1:39:40 PM
S100	11/19/2012 1:39:46 PM
UNC Mouse	11/19/2012 1:39:53 PM
UNC Rabbit	11/19/2012 1:40:00 PM
CD20	11/19/2012 11:16:14 AM
CD138	11/19/2012 11:27:36 AM
CD3	11/19/2012 11:16:21 AM
CD5 4C7	11/19/2012 11:16:28 AM
CD23 DAK-CD23	11/19/2012 11:16:35 AM
Kappa	11/19/2012 11:18:08 AM

Agilent Technologics



A Better Path[™]

🕜 DakoLink	A & 3 / 100	1000	1		- 0 X
System Administration Help					
New slides	Workflow	Instruments	Completed		
Completed sli	des				
Date completed: Any date	Case number	: Name:	Protocol:	Search Corinted Advar	Clear
To group by a column, drag its hea	ader here			View: *IHC default	•
Slide number Case number	Protocol	Primary antibody	Instrument Report printed	Completed Status	E
1	HercepTest	Her2 Ab	AS48	11/28/2012 2:02:53 PM Completed	
2	💽 Actin	Actin	AS48	11/28/2012 11:37:26 AN Completed	
3	HercepTest	Her2 NC	AS48	11/28/2012 2:02:53 PM Completed	
11	💽 Actin	Actin	AS48	11/28/2012 7:50:26 PM Completed	
···· 12	🔊 Alpha-1-Antitrypsin	Alpha-1-Antitrypsin	AS48	11/28/2012 7:50:26 PM Completed	
	💽 AMACR	AMACR	AS48	11/28/2012 7:50:26 PM Completed	
	BCL2	BCL2	AS48	11/28/2012 7:50:26 PM Completed	
- 15	🔊 BSAP	BSAP	AS48	11/28/2012 7:50:26 PM Completed	
- 16	🔊 Caldesmon	Caldesmon	AS48	11/28/2012 7:50:26 PM Completed	
	🔊 CD138	CD138	AS48	11/28/2012 7:50:26 PM Completed	
	🔊 CD1a	CD1a	AS48	11/28/2012 7:50:27 PM Completed	
== 24	 CD20	CD20	AS48	3/18/2013 11:13:34 AM Halted: E-Stop	
== 25	CD20	CD20	AS48	3/18/2013 11:13:34 AM Halted: E-Stop	
	CD20	CD20	AS48	3/18/2013 11:13:34 AM Halted: E-Stop	
	 CD20	CD20	AS48	3/18/2013 11:13:34 AM Halted: E-Stop	
	 CD20	CD20	AS48	3/18/2013 11:13:34 AM Halted: E-Stop	
3 0	CD10	CD10	AS48	3/18/2013 11:13:34 AM Halted: E-Stop	
		0010	1010		T
Print report	Refresh Search re	sults as of: 9/24/2014 3	:41:52 PM 105 slides		
Doyon, Dave [Administrator]	S: 127.0.0.1 C: 148.5.82.123	Site ID: 1000	DakoLink (3.0.34.987)		



🚺 DakoLink					
System Administration Help					
New slides Workflow	/ 🛒 In	struments	Completed		
Completed slides					
Date completed: Between ▼ 4/1/2014 ▼ - 4/22/2014 ▼	Case number:	Name:	Protocol:	Search	Clear Advanced
			Vie	ew: *HSC GRADING/SORT	ING 🔻
Silde number Case number Block Protoco + Instrument : HOWARD(AS3) (Slides = 881) + Instrument : LEONARD(AS2) (Slides = 948) + Instrument : SBH AS1 (Slides = 349) + Instrument : SBH AS2 (Slides = 384) + Instrument : SHELDON(AS1) (Slides = 1128) + Instrument : WL AS (Slides = 321)		Slide notes	Completed	▲ Target retrieval	
Print report Refresh	Search results a	as of: 4/22/2014 3:58:50 PM	4011 slides		
Doyon, Dave [Administrator] S: 127.0.0.1	C: 172.21.195.111	Site ID: 3000	DakoLink (3.0.34.987)		



Δ	Retter	Path™
/ \	DOLLOI	i aui

🚯 DakoLink				X
System Administration Help				-
	1			
New slides Wo	rkflow II	nstruments	Completed	
Completed slides				
Date completed:	Case number:	Name:	Protocol:	
Between 🕶 4/1/2014 💌 - 4/22/2014	•		CD20 Search Clear	
	🗌 My Slides		IHC Report not printed Advanced	
Instrument 🔺			View: *H5C GRADING/SORTING	
Slide number Case number Block	Protocol	Slide notes	Completed 🛆 Target retrieval	
Instrument : HOWARD(AS3) (Slides = 12)				
+ Instrument : LEONARD(AS2) (Slides = 33)				
+ Instrument : SBH AS1 (Slides = 8)				
+ Instrument : SBH AS2 (Slides = 7)				
+ Instrument : SHELDON(AS1) (Slides = 30)				
+ Instrument : WL AS (Slides = 5)				
Print report Ref	resh Search results	as of: 4/22/2014 4:02:27 PM	95 slides	
Doyon, Dave [Administrator] S: 127.0.0.1	C: 172.21.195.111	Site ID: 3000	DakoLink (3.0.34.987)	



Α	Better	Path [™]
	Delle	I aui

Fields:		Search criteria:			
Field	Δ 🔺	Field	Operator	Value	
Block					
Bottle serial					
Case number					
Completed					
Completed successfully	=				
Control slide					
Deparaffinization					
Drop zone					
Enzyme pretreatment					
External slide number					
First name					
Instrument					
Last name					
Lot					
Manual placed					
Middle name					
Drimany antibody	-				

Reports: Monitor usage



Slide Count by Protocol

Date from 1/1/2013

Date to 12/31/2013

Export to Excel

Primary antibody

Name	Description	Number of slides
AE1/3	FLEX Cytokeratin, MxH (AE1/AE3)	792
AFP	AFP, Rx	75
BCL-2	FLEX BCL2 Oncoprotein, MxH (124)	142
BCL-0		133
CALC		100
CALRE	Easier to:	220
CD 34		149
CD10	Find Root cause of staining issues	237
CD138	- Troublochocting with multi lob perspective	141
CD21 10	Troubleshooting with multi-lab perspective	2
CD21 30	Identify OC material	2
CD23		80
CD3	Understand tilization from 1 or multiple sites	2041
CD5 SP19		41
CD56	 Generate electronic reports for accreditations 	253
CD68		124
CDX-2	Manage Inter-lab QC proficiency	223
CHROM	- Deuteurs in deuth statistical data analysis	258
CK5/8	Perform in-depth statistical data analysis	309
CKIT		126
CYCD1	FLEX Cyclin D1, RxH (EP12)	14
Cyclin D1 10	Cyclin D1 10	2
Cyclin D1 20 Cyclin D1 20	Cyclin D1 20	2
CYT07	FLEX Cytokeratin 7 , MxH (OV-TL 12/30)	410
CYT20	FLEX Cytokeratin 20 , MxH (Ks20.8)	396
00.40	Total slides for	r group 25229



Wide Area Network (WAN) utilized applications at:



DIAGNOSTIC SERVICES

- LAN seat utilization
- Print on demand at Grace Hospital
- Centralized antibody preparation, optimization, and validations
- Seamless menu expansion and modifications
- Standardization of staining
- Multi-site contingency: Instruments, reagents
- Track and trace: Slides status
- Reports: Easy data mining, monitor usage
- Remote access customer support: Software upgrades, troubleshooting

Future connectivity solutions at DSM









Future connectivity solutions at DSM





- Addition of special stains in IHC database
- LIS interface
 - Eliminate double entry = save time
- Tracking solution



Suggestions for implementation



Decision Process

- Assess the level of tracking desired, define goals
- Check capabilities of current LIS system
- Does current LIS offer a tracking module?
- Does the module meet your goals?
 - If yes, implement
 - If no, seek other vendor middleware options
- Consider current internal IT capabilities before seeking outside help
 - Contract outside as needed



Suggestions for implementation



Must Do

- Clearly identify the purpose and goals of the project
- Involve your IT department from the start
- Learn about the limitations of your existing systems
- Outline capabilities needed from outside organization
- Look at solutions from multiple vendors
 - Don't assume that they all have the same capabilities
- Assign dedicated project manager
- Develop plan and implement in stages
 - Test system at key milestones
 - Catch and correct any problems early



Suggestions for implementation



Pitfalls to Avoid

- Lack of clearly defined goals
- Assuming all tracking systems are the same
- Taking on a project that is too large more than you need
- Not exploring compatibility of various systems, versions, etc.
- Involving IT department too late in the project
- Looking at connectivity as a side project
 - Lack of dedicated resource to manage



Summary of benefits - why connectivity?



- Eliminate manual data entry
- Reduce possibility of errors
- Maximize employee resources
- Track specimens and slides
- Manage workflow / workload
- Growth without increase in FTEs
- Better patient care!



Discussion











Anthony St-Onge – Strategic Business Manager Anthony.St-Onge@dako.com

Christopher Marino – Strategic Business Manager, Chris.marino@dako.com

