

Using Lean to Create High-Performance Workflow in a Histology R&D Laboratory

Dave Stanforth
Director Workflow Applications R&D
Dako Denmark A/S
Lab Quality Confab 2009

Agenda



- Why focus on workflow?
- Making the committment: building the Workflow Applications Lab
- Starting to get lean: how we use the lab
- Continuing the journey: where we're heading

Anatomic Pathology Laboratory Challenges

The Need to Improve Efficiencies

Cost and Productivity

- Increasing workload
- Labor shortages
- Budget constraints
- Labor-intensive work processes

Turn Around Time

- Pressure to reduce time to diagnosis
- Patient fast tracking (STAT requests)

(A) Dako DRAFT

Quality

- Specimen handling and tracking
- Data entry errors
- Staining quality and reproducibility



Workflow Learning Curve in AP

Market Segmentation Estimate



- Organize (Efficiency) 80% of labs
 - Remove non-value added and repetitive tasks.
 - Implement automation.
- Standardize (Quality) 10% of labs
 - Define and validate processes.
 - Document that validated processes are used.
 - Be alerted if there are deviations from the defined process.
- Optimize (Productivity) 10% of labs
 - Real time data collection to support LEAN lab metrics and KPI's.
 - Optimized and adaptive operations.
 - Proactive equipment maintenance and inventory management.

Lab Quality Confab 2009 DStanforth 4

Workflow Applications R&D




Dako's Concentrated Focus on Workflow

- Create a new department within the R&D function at the corporate headquarters.
- Expand our development focus from individual products to integrated solutions.
- Educate ourselves so that we can use lean tools and principles during product development, testing and validation.
- Use lean concepts to develop efficient standardized processes for Anatomic Pathology laboratories.
- Validate workflow solutions in a real-life Anatomic Pathology laboratory environment before launch.

Lab Quality Confab 2009 DStanforth

5

Anatomic Pathology Laboratory Typical Workflow Process Map Grossing Accessioning Sectioning Sectioning Staining Interpretation Lab Quality Contab 2009 DSanforth 6



Workflow Applications Laboratory AFT Floorplan – Not Lean!



Insert blueprints

Lab Quality Confab 2009 DStanforth

Homework: AP Lab Time Studies

Accessioning

1 Case with 2 containers



Observation	Time (s)	Cumulative (s)
Checking	60	60
Data entry	80	140
Labels on (6)	40	180
Time per container / Total time	90	180

Note: One case with two containers = 2 form labels + 4 container labels

Case entry to LIS (fast cases without clinic entry)


Observation	Time (s)	Cumulative (s)
Case 1	40	40
Case 2	28	68
Case 3	35	103
Case 4	43	146
Case 5	35	181
Case 6	44	225
Case 7	48	273
Case 8	24	297
Case 9	30	327
Case 10	30	357
Average per case / Total	35,7	357


After this 32s to print labels. One missing label was produced. That took 50s.


Lab Quality Confab 2009 DStanforth


Case entry into LIS for on-site cases (13 cases)


Observation	Time (s)	Cumulative (s)
Check and book clinics (13 cases)	137	137
Accession and print labels (13 cases)	355	492
Label 14 cases & 25 containers	184	676
Create 34 cassettes for the 25 containers	306	982
Avg check and clinics	10,5	
Avg accession and print	/27,3	
Avg labeling per case	13,1	
Avg labeling per container	7,4	
Afg cassette printing	9,0,	





























Workflow Applications Laborator PRAFT Dako

- Add H&E and Special Stains capabilities to enable the development of a total staining solution for complete cases.
- Revise floor plan to reduce unecessary travel and increase process visibility.
- Rebuild IT structure to enable efficient lab configuration changes.
- Support for distinct production and test environments.
- Establish performance benchmark processes to measure improvements.
- Optimize individual process step workflow.







Workflow Applications Laboratory Dako Improving Our Product Validation Process

- Validate new workflow solutions in a production environment by processing a realistic sample workload on a realistic schedule.
- Assure that all information and samples flow properly from start to finish and that we have well-defined processes for every step from accessioning to archiving.
- Create validation scenarios based on process mapping data collected at Anatomic Pathology labs.
- Map laboratory current state and validate proposed improvements utilizing new solutions and processes.

Lab Quality Confab 2009

Workflow Applications Laboratory Validation Every le

"Day in the Laboratory" Validation Example

- 07:00 Laboratory opens (6 pending cases)
- 07:30 Specimen Delivery (9 new cases)
- 7-11 Specimens from surgery (4 cases, one frozen)
- 11:00 Specimen Delivery (2 cases)
- 11-13 Specimens from surgery (2 cases)
- 13:30 Specimen Delivery (5 cases)
- 14-17 Specimens from surgery (1 case)
- 17:00 Laboratory closes

Lab Quality Confab 2009 DStanforth

25

